


Use of focus beads for imaging of cell-based assays


Quickly finding the correct focus level is one of the key challenges for automated, microscope based fluorescence imaging systems. In most applications DAPI-focussing on cell nuclei is used to determine the bottom of slides, plates or other carrier materials. However, this method can lead to errors when the cells detach from the bottom surface or are insufficiently stained.

PolyAn addresses this problem with our new 2 µm blue fluorescent Polymethylmethacrylate (PMMA) beads.


- 2 µm beads have a comparable size to cell nuclei.
- PMMA-grade is not cell-toxic and does not interfere with the cell behavior.
- Color encoded in the DAPI channel, transparent in all other channels.
- Well-defined fluorescence intensity. Problems caused by insufficient staining of the cell nuclei can thus be avoided.
- Easy handling: the beads can be easily added to the cell suspension. They quickly sink to the bottom and thus indicate the correct focus level.
- Application: standards in a range of cell assays, e.g. biofilm assays, adhesion assays or detection of bacteria.



In the examples illustrated below the read-out was done using the VideoScan HCU*, a fluorescence imaging system based on a fluorescence microscope**. The microscope focusses quickly on the focus beads independently whether bacteria adhere or not.



Focus beads (blue, DAPI channel) in combination with EPEC-bacteria on a GP2-coated plate. Only the surface of the bacteria has been stained with O26 *E.coli*-antibody sera and FITC-conjugated secondary antibody.



Focus beads (blue, DAPI channel) in combination with *E.coli*-bacteria on a GP2-coated plate. The bacteria have been stained using PI (DNA).

* Rödiger, S. et.al. Adv. Biochem. Eng. Biotechnol. 2013, 133, 35-74

** Schierack, P. et.al. Gut (BMJ Group), 2014, manuscript accepted, in press

High Performance Consumables

PolyAn is a nanotechnology company specialized in the modification of surfaces using Molecular Surface Engineering (MSE). Since 1996 PolyAn develops and manufactures consumables for multiplex diagnostic and LifeScience research.

Molecular Surface Engineering Services

PolyAn is able to equip almost any substrate with our reactive matrices for selective immobilization and antifouling surfaces for the reduction of cell adhesion and unspecific binding, respectively. As part of our Molecular Surface Engineering services, we offer functionalized consumable and substrate materials for OEM applications, which are tailored to specified customer requirements.

We are looking forward to your telephone orders and technical enquiries at our Customer Service and Technical Service Department Monday-Friday. Office hours for telephone enquiries are 9:00 AM to 6:00 PM (Central European Time). Please mention billing and shipping addresses, product-ID, quantity, your phone number or e-mail and name.

Distributors

**USA, Canada,
Mexico**

AutoMate Scientific, Inc. (USA)

Tel +1 (0)800 998.6283
+1 (0)510 845.6283
Fax +1 (0)510 280-3795
Email info@autom8.com

China

Beijing Thinkout Sci-Tech Co. Ltd.

Tel +86-13901326227
Email yuxb@pmbio.cn

**Singapore,
Malaysia**

Sciencewerke Pte. Ltd. (Singapore)

Tel +65 (0)6777 1045
Fax +65 (0)6777 3054
Email jason@sciencewerke.com

Korea

Kyongshin Scientific Co., Ltd.

Tel +82-2-576-6303
Fax + 82-2-576-6309
Email lee@kyongshin.co.kr

Taiwan

Bio-cando inc., Taiwan

Tel +886 (0)3 211-8079
Email bio-cando@promogene.com.tw

Israel

Moshe Stauber Biotec Applications

Tel +972-8-936 70 01
Fax +972-8-936 70 02
Email stauber1@zahav.net.il

Japan

CytoTechs, inc.

Tel +81-29-834-7788
Fax +81-29-834-7772
Email j.iijima@cytotechs.com

Italy

K.F. Technology Srl

Tel +39 06 454.34.179
Fax +39 069725.3131
Email fabrizio@kftechnology.it